






## American Procurement Services: Distributor of Noctic Medical Equipment



**American Procurement Service: Medical  
Equipment & Furniture Division**



**Please contact us for Pricing and product Information:  
Info@Americanprocurement.com or  
Sales@Americanprocurement.com**


Item	Pic. and Model	Name	Description	Specification
1	 APS-WC-01	Wheel Chair	spray coating soft seat with 16u	kid's crutch. Fabric seat, pneumatic tyre, fastern bell, with metal pedal, plastic pedal or aluminum pedal. Outside width 65, height 80, inner width 40, seat height 45. Steel material, available for 90kgs.

2	 <p>APS-WC-02</p>	<b>Wheel Chair</b>	aluminum soft seat with 24u	double brake, extra width front wheel, metal pedal, 90H*45W(In)*67L, seat height 45, aluminum material, available for 90kgs.
3	 <p>APS-WC-03</p>	<b>Wheel Chair</b>	spray coating toilet seat with 22u	fabric seat, fastern bell, 5' front wheel, 22' rare wheel, metal pedal, plastic holder, with toilet seat, 85H*45W(in)65W(out)*85L, steel material, available for 100kgs



4	 <p>APS-WC-04</p>	<p><b>Wheel Chair</b></p>	<p>aluminum folded holder soft seat with 16u</p>	<p>folded aluminum material, avialable for 90kgs</p>
5	 <p>APS-WC-05</p>	<p><b>Wheel Chair</b></p>	<p>aluminum toilet seat with 24u</p>	<p>aluminum material with toilet seat, double brakes. 90H*45W*95L, seat height 45, available for 100kgs,</p>

6	 <p data-bbox="353 639 495 667"><b>APS-WC-06</b></p>	<b>Wheel Chair</b>	plating with cushion open	fabric seat, fastern belt, plating steel material, extra width front wheel, 90H*45W*45H(seat), outside width 65, folded size: 92H*26W*90L, available for 100kgs
7	 <p data-bbox="353 1145 495 1173"><b>APS-WC-07</b></p>	<b>Wheel Chair</b>	plating soft seat	plastic handle, plating steel material, metal pedal with extra width front wheel. Folded size 90H*26W*90L, available for 90kgs


8	 <p>APS-WC-08</p>	Wheel Chair	plating toilet seat with non-air wheel	plating steel with toilet seat, non-slipe hand pulling wheeling, folded size 90H*26W*90L, available for 100kgs
9	 <p>APS-WC-09</p>	Wheel Chair	plating high back	ajusted lying degree, non-slipe hand pulling wheel with toilet seat, double brakes, fastern belt, folded size 95H*26W*95L, avialable for 100kgs

<p>10</p>	 <p>APS-WC-10</p>	<p><b>Wheel Chair</b></p>	<p>plating full lay down leather seat</p>	<p>leather seat, with toilet seat, adjusted laying degree, double brakes with non-slip pulling handle, available for 100kgs, folded size 95H*26W*93H</p>
<p>11</p>	 <p>APS-WC-11</p>	<p><b>Wheel Chair</b></p>	<p>aluminum with cushion opened, 24u</p>	<p>double brake, hard seat, folded size 90H*26W*95L, aluminum material, available for 100kgs</p>

12	 <p style="text-align: center;"><b>APS-WC-12</b></p>	<b>Wheel Chair</b>	spray with couison opened	black color, hard seat with extra cushion, plastic pedal, folded size 90H*26W*90L, steel material, available for 100kgs.
13	 <p style="text-align: center;"><b>APS-WC-13</b></p>	<b>Wheel Chair</b>	spray with toilet seat	silver color, extra width front wheel with toilet seat, non-slipe handle, metal pedel, steel material, available for 100kgs

14	 <p data-bbox="353 639 495 667"><b>APS-WC-14</b></p>	<b>Wheel Chair</b>	spray soft seat	silver color, steel material for 90kgs max. Folded size 90H*26W*90L
15	 <p data-bbox="353 1139 495 1166"><b>APS-WC-15</b></p>	<b>Wheel Chair</b>	spray soft seat with non-air wheel	black color with fabric seat, metal pulling wheel, folded size 90H*26W*90L, steel material for 90kgs max.



16	 <p data-bbox="353 639 495 667">APS-WC-16</p>	<b>Wheel Chair</b>	spray soft seat with non-air wheel	black color with fabric seat, metal pulling wheel, folded size 90H*26W*90L, steel material for 90kgs max.
----	--	--------------------	------------------------------------	---